



Variable information: LAS3B168 (questionnaire);
 LAS3B368 (technical measures)

LAS3B168

name	label	values
bmcat1	CAT1: cough	no answer, asked(-1)
bmcat2	CAT2: production of phlegm	no answer, asked(-1)
bmcat3	CAT3: chest tightness	no answer, asked(-1)
bmcat4	CAT4: breathlessness	no answer, asked(-1)
bmcat5	CAT5: activity limitation	no answer, asked(-1)
bmcat6	CAT6: confidence	no answer, asked(-1)
bmcat7	CAT7: sleep	no answer, asked(-1)
bmcat8	CAT8: energy	no answer, asked(-1)
bmcattot	CAT total score	no valid data(-2)
bmlftd	Long function: test done	no answer, asked(-1), no(-1), yes(2)
bmlftrn	Long function: reason test not done	LF test done(-2), device not working or parts missing(1), physical reasons(2), R does not want to perform spirometry(3), limited time(4), no medication(5), environment not suitable(6), unknown(7)

LAS3B368

name	label	values
FEV1PREDQuanjerGLIreference	FEV1PRED (Quanjer, GLI reference)	none
FEV1liters	FEV1 liters	none
FVCliters	FVC liters	none
Tiff		none
FEV1LLN	FEV1 LLN	none
FVCLLN	FVC LLN	none
FEV1FVCLLN	FEV1/FVC LLN	none
Tiff_LL_NHANES		none
COPD_LL_NHANES		none
Tiff<70	Tiff<70	none
StageIV		none
StageIII		none
StageII		none
StageI		none
SpirometryError	Spirometry Error	none
Insufficienttotalerror	Insufficient total error	none
FEV1error	FEV1 error	none
FVCerror	FVC error	none
Expirationerror		none
Practicesession		none