



Variable information: LASAB172

Bioelectric impedance is not measured in the whole sample, only in interviews taken in 1993 (N=247).

name	label	values
bmresist	Bioelectric impedance: resistance	not in sample(-3), no valid measurement(-1)
bmreacta	Bioelectric impedance: reactance	not in sample(-3), no valid measurement(-1)
bmrefuse	Respondent refused	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), reason not applicable(1), reason applicable(2)
bmno plac	Not an adequate place to lie down	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), reason not applicable(1), reason applicable(2)
bmphysic	Respondent not physically able to undergo the measurement	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), reason not applicable(1), reason applicable(2)
bmequipm	Equipment malfunction	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), reason not applicable(1), reason applicable(2)
bmedema	R characteristics: edema on leg	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), no(1), yes(2)
bm cold	R characteristics: cold (shivering)	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), no(1), yes(2)
bmdrink	R drank before measurement	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1), coffee(1), tea(2), other(3), did not drink(4)
bmnumber	R number of drinks	not in sample(-3), no answer, skipped(-2), no answer, coding error(-1),