

Variable information: LAS3B168 (questionnaire);  
LAS3B368 (technical measures)

**LAS3B168**

name	label	values
bmcat1	CAT1: cough	no answer, asked(-1)
bmcat2	CAT2: production of phlegm	no answer, asked(-1)
bmcat3	CAT3: chest tightness	no answer, asked(-1)
bmcat4	CAT4: breathlessness	no answer, asked(-1)
bmcat5	CAT5: activity limitation	no answer, asked(-1)
bmcat6	CAT6: confidence	no answer, asked(-1)
bmcat7	CAT7: sleep	no answer, asked(-1)
bmcat8	CAT8: energy	no answer, asked(-1)
bm cattot	CAT total score	no valid data(-2)
bm lftd	Long function: test done	no answer, asked(-1), no(-1), yes(2)
bm lftfn	Long function: reason test not done	LF test done(-2), device not working or parts missing(1), physical reasons(2), R does not want to perform spirometry(3), limited time(4), no medication(5), environment not suitable(6), unknown(7)

**LAS3B368**

name	label	values
FEV1PREDQuanjerGLIreference	FEV1PRED (Quanjer, GLI reference)	none
FEV1liters	FEV1 liters	none
FVCliters	FVC liters	none
Tiff		none
FEV1LLN	FEV1 LLN	none
FVCLLN	FVC LLN	none
FEV1FVCLLN	FEV1/FVC LLN	none
Tiff_LL_NHANES		none
COPD_LL_NHANES		none
Tiff<70	Tiff<70	none
StageIV		none
StageIII		none
StageII		none
StageI		none
SpirometryError	Spirometry Error	none
Insufficienttotalerror	Insufficient total error	none
FEV1error	FEV1 error	none
FVCerror	FVC error	none
Expirationerror		none
Practicesession		none