

Variable information: LASAC185 / LASAD185 / LASAE185 / LASAF185

NB.: In wave C, ankle, knee and pulse width (see questionnaire LASAC161) were processed in LASAC185.

LASAC185 / LASAD185

name	label	values
c/dmfeetinsp	Inspection feet: whole	-5: interview terminated (<i>only in C</i>) -1: no valid data 1: yes 2: no
c/dmfeetinsp1	Inspection feet: sores	-2: na, see C/DMFEETINSP
c/dmfeetinsp2	Inspection feet: plaster	-1: no valid data / na, asked
c/dmfeetinsp3	Inspection feet: corns	0: not mentioned
c/dmfeetinsp4	Inspection feet: skewed toes	1: mentioned
c/dmfeetinsp5	Inspection feet: amputation toes	
c/dmfeetinsp6	Inspection feet: amputation part of foot	
c/dmfeetinsp7	Inspection feet: amputation feet	
c/dmfeetinsp8	Inspection feet: other	
c/dmfeetinspo	Inspection feet: other specification	-3: na, wrong skip -2: na, see C/DMFEETINSP8 -1: no valid data / na, asked 0: to be coded
c/dmfeetpllts	Inspection feet: pull up toes	-5: interview terminated (<i>only in C</i>)
c/dmfeetsmf	Inspection feet: smack feet	-1: no valid data / na, asked 1: yes 2: no
Only in C:		
cmankle1	Ankle (1): width in mm	-5: na, interview terminated
cmankle2	Ankle (2): width in mm	-1: no valid data
cmknee1	Knee (1): width in mm	
cmknee2	Knee (2): width in mm	
cmpulse1	Pulse (1): width in mm	
cmpulse2	Pulse (2): width in mm	

LASAE185 / LASAF185

name	label	values
e/fmfeetpllts	Inspection feet: pull up toes	-1: na, asked 1: yes 2: no